



Large Orifice (Wide Bore) Tips



- Larger opening prevents cell fragmentation due to shearing.
- Ideal for pipetting of macromolecules and viscous samples.
- Inside diameter at tip end is 0.06 in (1.6 mm).
- Free of detectable RNase, DNase, DNA, and PCR inhibitors. Sterile versions are also pyrogen-free.

Item #	Volume	Description	Color	Packaging
4290-00	200 µL	Large Orifice Tips, Bulk	Natural	1000 tips/unit, 10 units/case
4297-00	200 µL	Large Orifice Tips, Racked	Natural	10 racks/unit, 5 units/case
4297-S0	200 µL	Large Orifice Tips, Racked, Sterile	Natural	10 racks/unit, 5 units/case
4390-00	1 mL	Large Orifice Tips, Bulk	Natural	1000 tips/unit, 10 units/case
4397-00	1 mL	Large Orifice Tips, Racked	Natural	10 racks/unit, 5 units/case
4397-S0	1 mL	Large Orifice Tips, Racked, Sterile	Natural	10 racks/unit, 5 units/case

Interpath Services
Leading the way

Available from: Interpath Services
www.interpath.com.au
Ph: (03) 9457 6277



Made in USA