





Interpath Services[®]
Leading the way



Aerosol Barrier Tips



Aerosol Barrier Tips







SIZE	PART NUMBER	COMPATIBLE WITH						
		Socorex	Biohit	Brand Transferpette	Eppendorf	Finnpipette	Gilson Pipetman	Nichiryo/Oxford
10	24100 Low Retention Barrier Tip	Acura 8/12 Channel 0.5 - 10µl Acura 0.5 - 10µl	Electronic 0.2 - 10µl Proline 0.5 - 10µl	0.1 - 10µl	Electronic 0.5 - 10µl Reference 0.1 - 10µl 0.5 - 10µl Research 0.5 - 10µl		P2 P10	Benchmate 0.5 - 10µl
	Tip length: 31.25mm							
10 F	24200 Low Retention Barrier Tip	Acura (Y Series) 1 - 10µl	Proline 0.5 - 10µl			Single 0.5 - 20µl 5 - 40µl		
	Tip length: 50mm							
10 XL	24300 Low Retention Barrier Tip	Acura 0.1 - 2µl 0.5 - 10µl Acura 8/12 Channel 0.5 - 10µl	Electronic 0.2 - 10µl Proline 0.5 - 10µl	0.5 - 10µl	Reference 0.1 - 2.5µl 0.5 - 10µl Research 0.1 - 2.5µl 0.5 - 10µl	0.5 - 10µl	P2 P10	Benchmate 0.1 - 2µl 0.5 - 10µl
	Tip length: 45.84mm							
10 E	24400 Low Retention Barrier Tip		Proline 0.5 - 10µl	0.1 - 10µl	Reference 0.1 - 10µl Research 0.1 - 2.5µl 0.5 - 10µl	Single 0.5 - 20µl	P2 P10	Benchmate 0.5 - 10µl
	Tip length: 45mm							

Tip photographs not actual size.

All tips up to 1000µl are packed in 10 trays of 96 tips. The 1000µl and 1250µl tips are packed in 8 trays of 96 tips.



Certified
RNase/DNase/Pyrogen
FREE

SIZE	PART NUMBER		COMPATIBLE WITH						
			Socorex	Biohit	Brand Transferpette	Eppendorf	Finnpipette	Gilson Pipetman	Nichiryo/ Oxford
20	24500	Low Retention Barrier Tip	Acura 2 - 20µl Acura 8/12 Channel 5 - 50µl	Electronic 0.2 - 10µl 5 - 100µl Proline 5 - 50µl 50 - 200µl	5 - 50µl 20 - 200µl	Reference 2 - 20µl Research 2 - 20µl	Single 5 - 40µl Multi- channel 5 - 50µl	P20	Benchmate 2 - 20µl
		Tip length: 49.63mm							
100	24600	Low Retention Barrier Tip	Acura 10 - 100µl 20 - 200µl Acura 8/12 Channel 20 - 200µl 40 - 350µl	Electronic 0.2 - 10µl 5 - 100µl Proline 5 - 50µl 50 - 200µl	5 - 50µl 50 - 100µl 20 - 200µl	Electronic 5 - 100µl 2 - 200µl Reference 10 - 100µl 20 - 200µl Research 10 - 100µl 20 - 200µl	Single 40 - 200µl	P100	Benchmate 10 - 100µl 20 - 200µl
		Tip length: 59.54mm							
200	24700	Low Retention Barrier Tip	Acura 20 - 200µl Acura 8/12 Channel 20 - 200µl 40 - 350µl	Electronic 5 - 100µl Proline 5 - 50µl 50 - 200µl	10 - 100µl 20 - 200µl	Electronic 5 - 100µl 2 - 200µl Reference 10 - 100µl 20 - 200µl Research 10 - 100µl 20 - 200µl	Single 40 - 200µl Multi- channel 50 - 300µl	P200	Benchmate 10 - 100µl 20 - 200µl
		Tip length: 59.54mm							
300	24750	Low Retention Barrier Tip	Acura 20 - 200µl	Electronic 5 - 100µl 20 - 200µl Proline 10 - 100µl 20 - 200µl	20 - 200µl 15 - 300µl	Reference 50 - 200µl Research 20 - 200µl 30 - 300µl	Single 5 - 50µl 20 - 200µl 30 - 300µl Multi- channel 5 - 50µl 30 - 300µl	P200	Benchmate 20 - 200µl 50 - 300µl
		Tip length: 59.44mm							
1000	24800	Low Retention Barrier Tip	Acura 100 - 1000µl	Proline 200 - 1000µl	100 - 1000µl	Electronic 50 - 1000µl Reference 100 - 1000µl Research 100 - 1000µl	Single 100 - 1000µl 200 - 1000µl	P1000	Benchmate 100 - 1000µl
		Tip length: 83.94mm							
1250	24900	Low Retention Barrier Tip	Acura 100 - 1000µl	Proline 200 - 1000µl	100 - 1000µl	Reference 100 - 1000µl Research 100 - 1000µl	Single 100 - 1000µl 200 - 1000µl Multi- channel 100 - 1200µl	P1000	Benchmate 100 - 1000µl
		Tip length: 98.30mm							

Aerosol Barrier Tips

Interpath's Aerosol Barrier Tips are produced from the highest grade polypropylene polymer together with a unique virgin HDPE (High Density Polyethylene) barrier.

FEATURES	BENEFITS
<ul style="list-style-type: none"> ▼ Universal soft fit style 	<ul style="list-style-type: none"> ▼ Easy to acquire and eject tips ▼ Will suit a much wider range of pipettes than other brands
<ul style="list-style-type: none"> ▼ Much larger headspace 	<ul style="list-style-type: none"> ▼ Eliminates over pipetting which can result in possible contamination
<ul style="list-style-type: none"> ▼ Reinforced see-through hinged lid tip rack with security clasp 	<ul style="list-style-type: none"> ▼ Clear hinged lid racks allow for instant tip inspection, also convenience and ease of use ▼ Enhanced tip support when using multichannel pipettes ▼ Eliminates tip movement during use and transport
<ul style="list-style-type: none"> ▼ Low Retention tips 	<ul style="list-style-type: none"> ▼ Slippery tip surface to reduce sample retention and improve volume accuracy
<ul style="list-style-type: none"> ▼ 0.2µl barrier pore size 	<ul style="list-style-type: none"> ▼ Smaller pore size guarantees maximum possible aerosol protection. No additives or gums are used (as found in some brands – which can cause possible contamination and inhibit the PCR process). <i>Refer to important note below</i>
<ul style="list-style-type: none"> ▼ Highest manufacturing standards 	<ul style="list-style-type: none"> ▼ Tips are manufactured in accordance with ISO 9001 Standard, are FDA approved and certified RNase, DNase & Pyrogen Free
<ul style="list-style-type: none"> ▼ Lot numbers printed on individual racks 	<ul style="list-style-type: none"> ▼ Improved tip traceability



IMPORTANT NOTE:

The barrier of the Aerosol Barrier Tips have the smallest pore size (0.22µm) to ensure the most minute aerosol particles are trapped, thus preventing cross contamination. Unlike some brands, Interpath's Aerosol Barrier Tips do not contain a cellulose gum additive. Other brand tips using this gum have a larger (0.45µm) pore size and therefore offer greater risk of contamination as aerosol particles can pass through this pore size. Interpath's Aerosol Barrier Tips **contain no PCR inhibitors** and are tested using both human and bacterial templates. They are quality tested and contain less than 5×10^{-11} µg of human or bacterial DNA, less than 3.4×10^{-11} µg RNase and less than 1.7×10^{-11} µg DNase.

Contact Interpath Services for additional information

Unit 1 46 Sheehan Road
West Heidelberg VIC 3081
Ph: 03 9457 6277
Fax: 03 9458 4010

Freecall: 1800 62 63 69
Freefax: 1800 65 82 10
E-Mail: sales@interpath.com.au

Interpath Services[®]
Leading the way